


CRIO

CRYO RANGE

CRYOGENIC HANDLING

- Cryogenic handling in industry, pharmaceuticals, chemical sector, laboratories



AET N°0072/011/162/07/19/0268 (IFTH LYON)

Compliance with immersion in liquid nitrogen (IFTH MR 019 method)

STANDARDS (DETAILS ON WEB SITE)
www.rostaing.com/normes-secureite

SIZES



COMPOSITION

Palm and back Water-repellent cow leather
Cuff Cow split leather
Membrane Porelle®
Lining Wool/cotton fleece and ISO AIR felt
Total length 40 cm

CARE



PACKAGING

Bag of 5 pairs - Box of 20 pairs

PRO•LINE:
Our selection



Resistance to cryogenic cold

Composition resistant to splashes of liquid nitrogen at -196° C



Optimised service life

Reinforcement on the wear zone



Waterproof

Water-repellent leather does not harden once dry



Cuff length 15 cm



Adjustable strap

CRYO

